




LMO4 Antibody

Product Code	CSB-PA118935
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61968
Immunogen	Synthesized peptide derived from C-terminal of Human LMO4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total LMO4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	LIM domain only 4; LIM domain transcription factor LMO4; LMO-4; breast tumor autoantigen;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LMO4

Image



Western blot analysis of extracts from 293 cells, using LMO4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.