




IRX3 Antibody

Product Code	CSB-PA090661
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78415
Immunogen	Synthesized peptide derived from internal of Human IRX3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total IRX3 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	iroquois homeobox 3; IRX-1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IRX3

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using IRX3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.