




ST5 Antibody

Product Code	CSB-PA195797
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78524
Immunogen	Synthesized peptide derived from internal of Human ST5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ST5 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Suppression of tumorigenicity 5; HeLa tumor suppression 1; DENN domain-containing protein 2B; ST5; DENND2B
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ST5

Image



Western blot analysis of extracts from COLO205 cells, using ST5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.