




ADARB1 Antibody

Product Code	CSB-PA949783
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78563
Immunogen	Synthesized peptide derived from internal of Human ADARB1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ADARB1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 3.5.-.-; dsRNA adenosine deaminase; RNA-editing deaminase 1; RNA-editing enzyme 1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADARB1

Image



Western blot analysis of extracts from HepG2 cells, using ADARB1 antibody.