





RBP5 Antibody

Product Code	CSB-PA907983
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82980
Immunogen	Synthesized peptide derived from internal of Human CRBP III.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CRBP III protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Retinol-binding protein III; cellular; CRBP-III; HRBPiso; RBP5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RBP5

Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using CRBP III antibody.



Immunofluorescence analysis of LOVO cells, using CRBP III antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.