





TFAM Antibody

Product Code	CSB-PA003208
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00059
Immunogen	Synthesized peptide derived from internal of Human TFAM.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TFAM protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	RTEF1; TCF13L1; TEA domain family member 4; TEAD-4; TEAD4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TFAM

Image



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue, using TFAM antibody.



Immunofluorescence analysis of HepG2 cells, using TFAM antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.