




AKAP17A Antibody

Product Code	CSB-PA024951
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02040
Immunogen	Synthesized peptide derived from C-terminal of Human SFRS17A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SFRS17A protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	A-kinase anchor protein 17A; 721P; B-lymphocyte antigen; Protein XE7; Protein kinase A-anchoring protein 17A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP17A

Image



Western blot analysis of extracts from NIH-3T3 cells, using SFRS17A antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.