




SP3/SP4 Antibody

Product Code	CSB-PA081786
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02447
Immunogen	Synthesized peptide derived from C-terminal of Human SP3/4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total SP3/4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	LAG1 longevity assurance; SP260; TMSG1; Tumor metastasis-suppressor gene 1 protein;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SP3/SP4

Image



Western blot analysis of extracts from Jurkat cells, using SP3/4 antibody.