




CSTF1 Antibody

Product Code	CSB-PA255640
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05048
Immunogen	Synthesized peptide derived from N-terminal of Human CSTF1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CSTF1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CF-1 50 kDa subunit; CSTF 50 kDa subunit; CSTF1; Cleavage stimulation factor; 50 kDa subunit
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CSTF1

Image



Western blot analysis of extracts from 293 cells and HUVEC cells, using CSTF1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.