




CHAF1A Antibody

Product Code	CSB-PA943110
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13111
Immunogen	Synthesized peptide derived from internal of Human CAF1A.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total CAF1A protein.
Tested Applications	ELISA, WB, IHC, IF; WB: 1:500-1:3000, IHC: 1:50-1:200, IF: 1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Chromatin assembly factor 1 subunit A; CAF-1 subunit A; Chromatin assembly factor I p150 subunit; CAF-I 150 kDa subunit; CAF-I p150
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHAF1A

Image



Western blot analysis of extracts from K562 cells and A549 cells, using CAF1A antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.