




PRKG2 Antibody

Product Code	CSB-PA826106
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13237
Immunogen	Synthesized peptide derived from internal of Human CGK 2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CGK 2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	cGMP-dependent protein kinase 2; EC 2.7.11.12; CGK 2; Cgkii; Type II cGMP-dependent protein kinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKG2

Image



Western blot analysis of extracts from HUVEC cells, using CGK 2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.