




SMAD4 Antibody

Product Code	CSB-PA825890
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13485
Immunogen	Synthesized peptide derived from internal of Human Smad4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Smad4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	SMAD 4; Mothers against decapentaplegic homolog 4; Mothers against DPP homolog 4; Deletion target in pancreatic carcinoma 4; hSMAD4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMAD4

Image



Western blot analysis of extracts from HeLa cells, COS-7 cells, HepG2 cells and HUVEC cells, using Smad4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.