





TRA2A Antibody

Product Code	CSB-PA082967
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13595
Immunogen	Synthesized peptide derived from internal of Human TRA-2 α .
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TRA-2 α protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transformer-2 protein homolog; TRA-2 alpha;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRA2A


Image



Western blot analysis of extracts from HUVEC cells, using TRA-2 α antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TRA-2 α antibody.



Immunofluorescence analysis of HeLa cells, using TRA-2 α antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.