




TNFAIP1 Antibody

Product Code	CSB-PA094566
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13829
Immunogen	Synthesized peptide derived from internal of Human TNAP1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TNAP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B12; B61; EDP1; MGC2317; TNAP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFAIP1

Image



Western blot analysis of extracts from HUVEC cells treated with PMA (125ng/ml, 30mins), and HUVEC cells treated with serum (20%, 30mins), using TNAP1 antibody.