





CBFB Antibody

Product Code	CSB-PA107678
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13951
Immunogen	Synthesized peptide derived from internal of Human CBF β.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CBF β protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Core-binding factor subunit beta; CBF-beta; Polyomavirus enhancer-binding protein 2 beta subunit; PEBP2-beta; PEA2-beta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CBFB

Image



Western blot analysis of extracts from HeLa cells and HUVEC cells, using CBF β antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using CBF β antibody.