




E2F2 Antibody

Product Code	CSB-PA786782
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14209
Immunogen	Synthesized peptide derived from internal of Human E2F2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total E2F2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription factor E2F2; E2F-2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	E2F2

Image



Western blot analysis of extracts from HepG2 cells, K562 cells and COLO205 cells, using E2F2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.