




GCKR Antibody

Product Code	CSB-PA276533
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14397
Immunogen	Synthesized peptide derived from internal of Human Glucokinase Regulator.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Glucokinase Regulator protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Glucokinase regulatory protein; Glucokinase regulator;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GCKR

Image



Western blot analysis of extracts from mouse heart cells and rat heart cells, using Glucokinase Regulator antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.