




GK3P Antibody

Product Code	CSB-PA156631
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14409
Immunogen	Synthesized peptide derived from internal of Human GLPK3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GLPK3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Glycerol kinase; testis specific 1; EC 2.7.1.30; ATP:glycerol 3-phosphotransferase; Glycerokinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GK3P

Image



Western blot analysis of extracts from 293 cells, using GLPK3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.