




NR1I3 Antibody

Product Code	CSB-PA586225
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14994
Immunogen	Synthesized peptide derived from internal of Human NR1I3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NR1I3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CAR; Constitutive androstane receptor; Orphan nuclear receptor MB67; Orphan nuclear receptor NR1I3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NR1I3

Image



Western blot analysis of extracts from Jurkat cells, treated with serum (20%, 15mins), using NR1I3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.