




IRF4 Antibody

Product Code	CSB-PA250376
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15306
Immunogen	Synthesized peptide derived from internal of Human IRF4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total IRF4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	interferon regulatory factor 4; LSIRF; lymphocyte-specific interferon regulatory factor; multiple myeloma oncogene 1; MUM1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IRF4

Image



Western blot analysis of extracts from HeLa cells and Jurkat cells, using IRF4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.