




ANKRD1 Antibody

Product Code	CSB-PA297327
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15327
Immunogen	Synthesized peptide derived from internal of Human ANKRD1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ANKRD1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ALRP; C-193; C193; CARP; CVARP
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANKRD1

Image



Western blot analysis of extracts from COLO205 cells, using ANKRD1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.