




EPHA7 Antibody

Product Code	CSB-PA235595
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15375
Immunogen	Synthesized peptide derived from internal of Human EPHA7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total EPHA7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ephrin type-A receptor 7 [Precursor]; EC 2.7.10.1; Tyrosine-protein kinase receptor EHK-3; EPH homology kinase 3; Receptor protein-tyrosine kinase HEK11
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPHA7

Image



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using EPHA7 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.