




TAF13 Antibody

Product Code	CSB-PA775954
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15543
Immunogen	Synthesized peptide derived from internal of Human TAF13.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total TAF13 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription initiation factor TFIID subunit 13; Transcription initiation factor TFIID 18 kDa subunit; TAF(II)18; TAFII-18; TAFII18
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TAF13

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using TAF13 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.