




TAF1A Antibody

Product Code	CSB-PA949229
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15573
Immunogen	Synthesized peptide derived from internal of Human TAF1A.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total TAF1A protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TATA box-binding protein-associated factor 1A; TBP-associated factor 1A; TBP-associated factor RNA polymerase I 48 kDa; TAFI48; Transcription factor SL1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TAF1A

Image



Western blot analysis of extracts from HeLa cells, using TAF1A antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.