





ZNF174 Antibody

Product Code	CSB-PA546808
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15697
Immunogen	Synthesized peptide derived from internal of Human ZNF174.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZNF174 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein 174; AW-1; Zinc finger and SCAN domain-containing protein 8; ZNF174; ZSCAN8
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF174

Image



Western blot analysis of extracts from HeLa cells, using ZNF174 antibody.



Immunofluorescence analysis of HUVEC cells, using ZNF174 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.