





MAPKAPK3 Antibody

Product Code	CSB-PA247267
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16644
Immunogen	Synthesized peptide derived from internal of Human MAPK3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MAPK3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	MAP kinase-activated protein kinase 3; EC 2.7.11.1; MAPK-activated protein kinase 3; MAPKAP kinase 3; MAPKAPK-3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAPKAPK3

Image



Western blot analysis of extracts from COLO205 cells, using MAPK3 antibody.



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using MAPK3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.