




NFIL3 Antibody

Product Code	CSB-PA093754
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16649
Immunogen	Synthesized peptide derived from internal of Human NFIL3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total NFIL3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Nuclear factor interleukin-3-regulated protein; E4 promoter binding-protein 4; Interleukin-3 promoter transcriptional activator; Transcriptional activator NF-IL3A; Interleukin-3 binding-protein 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFIL3

Image



Western blot analysis of extracts from 293 cells and K562 cells, using NFIL3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.