





MAPK6 Antibody

Product Code	CSB-PA073951
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16659
Immunogen	Synthesized peptide derived from internal of Human p97 MAPK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total p97 MAPK protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 6; EC 2.7.11.24; Extracellular signal-regulated kinase 3; ERK-3; MAP kinase isoform p97
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK6

Image



Western blot analysis of extracts from K562 cells, using p97 MAPK antibody.



Immunofluorescence analysis of MCF-7 cells, using p97 MAPK antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.