




DAK Antibody

Product Code	CSB-PA184114
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3LXA3
Immunogen	Synthesized peptide derived from internal of Human DAK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total DAK protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Dihydroxyacetone kinase; EC 2.7.1.29; Glycerone kinase; DHA kinase; DAK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TKFC

Image



Western blot analysis of extracts from K562 cells and A549 cells, using DAK antibody.