




ADCK5 Antibody

Product Code	CSB-PA186715
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3MIX3
Immunogen	Synthesized peptide derived from internal of Human ADCK5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ADCK5 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Uncharacterized aarF domain-containing protein kinase 5; EC 2.7.11.-; ADCK5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADCK5

Image



Western blot analysis of extracts from HeLa cells, 293 cells and HUVEC cells, using ADCK5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.