




TP53I3 Antibody

Product Code	CSB-PA563966
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53FA7
Immunogen	Synthesized peptide derived from internal of Human QORX.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total QORX protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Shugoshin-like 1; hSgo1; Serologically defined breast cancer antigen NY-BR-85; SGO1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TP53I3

Image



Western blot analysis of extracts from 293 cells, using QORX antibody.