




DGKK Antibody

Product Code	CSB-PA099899
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5KSL6
Immunogen	Synthesized peptide derived from C-terminal of Human DGKK.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DGKK protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.1.107; Diglyceride kinase kappa; DGK-kappa; DAG kinase kappa; 142 kDa diacylglycerol kinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DGKK

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using DGKK antibody.



Immunofluorescence analysis of COS-7 cells,
using DGKK antibody.