





CDKL4 Antibody

Product Code	CSB-PA778730
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5MAI5
Immunogen	Synthesized peptide derived from C-terminal of Human CDKL4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CDKL4 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.22; Cyclin-dependent kinase-like 4; CDKL4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDKL4

Image



Western blot analysis of extracts from HepG2 cells and COLO205 cells, using CDKL4 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using CDKL4 antibody.



Immunofluorescence analysis of HeLa cells, using CDKL4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.