




TTC23 Antibody

Product Code	CSB-PA087958
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5W5X9
Immunogen	Synthesized peptide derived from internal of Human TTC23.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TTC23 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TPR repeat protein 23; Tetratricopeptide repeat protein 23; Cervical cancer proto-oncogene 8 protein; HCC-8; HCC8
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TTC23

Image



Western blot analysis of extracts from RAW264.7 cells, using TTC23 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.