




SFXN4 Antibody

Product Code	CSB-PA253040
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P4A7
Immunogen	Synthesized peptide derived from N-terminal of Human SFXN4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SFXN4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Breast cancer resistance marker 1; Sideroflexin-4; BCRM1; SFXN4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SFXN4

Image



Western blot analysis of extracts from HUVEC cells, using SFXN4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.