




SFXN4 Antibody

Product Code	CSB-PA253040
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P4A7
Immunogen	Synthesized peptide derived from N-terminal of Human SFXN4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SFXN4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Breast cancer resistance marker 1; Sideroflexin-4; BCRM1; SFXN4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SFXN4

Image



Western blot analysis of extracts from HUVEC cells, using SFXN4 antibody.