




TNFAIP8L2 Antibody

Product Code	CSB-PA285296
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P589
Immunogen	Synthesized peptide derived from internal of Human TNFAIP8L2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TNFAIP8L2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor; alpha-induced protein 8-like protein 2; TNFAIP8-like protein 2; TIPE2; Inflammation factor protein 20
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFAIP8L2

Image



Western blot analysis of extracts from HuvEc cells, using TNFAIP8L2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.