





ARMCX2 Antibody

Product Code	CSB-PA591675
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L311
Immunogen	Synthesized peptide derived from internal of Human ARMX2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ARMX2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Armadillo repeat-containing X-linked protein 2; Protein ALEX2; ARM protein lost in epithelial cancers on chromosome X 2; ARM CX2; ALEX2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARM CX2

Image



Western blot analysis of extracts from 293 cells, using ARMX2 antibody.



Immunofluorescence analysis of HepG2 cells,
using ARMX2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.