




MAEA Antibody

Product Code	CSB-PA203871
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L5Y9
Immunogen	Synthesized peptide derived from internal of Human MAEA.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MAEA protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Macrophage erythroblast attacher; Erythroblast macrophage protein; Cell proliferation-inducing gene 5 protein; Lung cancer-related protein 10; HLC-10
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAEA

Image



Western blot analysis of extracts from COLO205 cells, using MAEA antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.