




CTU1 Antibody

Product Code	CSB-PA992421
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7A3
Immunogen	Synthesized peptide derived from C-terminal of Human ATBP3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ATBP3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cancer-associated gene protein; ATP-binding domain-containing protein 3; ATPBD3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTU1

Image



Western blot analysis of extracts from LOVO cells, using ATBP3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.