




ADCK1 Antibody

Product Code	CSB-PA562737
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TW2
Immunogen	Synthesized peptide derived from internal of Human ADCK1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ADCK1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Uncharacterized aarF domain-containing protein kinase 1 [Precursor]; EC 2.7.11.-; ADCK1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADCK1

Image



Western blot analysis of extracts from A549 cells and RAW264.7 cells, using ADCK1 antibody.