




PLXDC1 Antibody

Product Code	CSB-PA548726
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUUK5
Immunogen	Synthesized peptide derived from Internal of Human PLXDC1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PLXDC1 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Plexin domain-containing protein 1; Tumor endothelial marker 7; Tumor endothelial marker 3; TEM3; TEM7
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLXDC1

Image



Immunohistochemical analysis of paraffin-embedded human brain tissue using PLXDC1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.