





TUSC5 Antibody

Product Code	CSB-PA062014
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXB3
Immunogen	Synthesized peptide derived from internal of Human TUSC5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TUSC5 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor suppressor candidate 5; Protein located at seventeen-p-thirteen point three 1; TUSC5; LOST1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUSC5

Image



Western blot analysis of extracts from LOVO cells, using TUSC5 antibody.



Immunofluorescence analysis of HeLa cells, using TUSC5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.