





HIPK4 Antibody

Product Code	CSB-PA999254
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NE63
Immunogen	Synthesized peptide derived from internal of Human HIPK4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HIPK4 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Homeodomain-interacting protein kinase 4; EC 2.7.11.1; HIPK4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HIPK4

Image



Western blot analysis of extracts from HUVEC cells, using HIPK4 antibody.



Immunohistochemistry analysis of paraffin-embedded human testis tissue using HIPK4 antibody.



Immunofluorescence analysis of A549 cells, using HIPK4 antibody.