




ZNF397 Antibody

Product Code	CSB-PA697767
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NF99
Immunogen	Synthesized peptide derived from N-terminal of Human ZNF397.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ZNF397 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Zinc finger protein 397; Zinc finger and SCAN domain-containing protein 15; Zinc finger protein 47; ZNF397; ZNF47
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF397

Image



Western blot analysis of extracts from HeLa cells and HUVEC cells, using ZNF397 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.