




ANKK1 Antibody

Product Code	CSB-PA995860
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFD2
Immunogen	Synthesized peptide derived from internal of Human ANKK1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ANKK1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ankyrin repeat and protein kinase domain-containing protein 1; EC 2.7.11.1; Protein kinase PKK2; X-kinase; Sugen kinase 288
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANKK1

Image



Western blot analysis of extracts from HeLa cells, using ANKK1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.