





ARID1B Antibody

Product Code	CSB-PA210417
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFD5
Immunogen	Synthesized peptide derived from interal of Human BAF250B.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total BAF250B protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ARID domain-containing protein 1B; Osa homolog 2; hOsa2; p250R; BRG1-binding protein hELD/OSA1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARID1B

Image



Western blot analysis of extracts from LOVO cells, using BAF250B antibody.



Immunofluorescence analysis of HUVEC cells,
using BAF250B antibody.