




ADCK3 Antibody

Product Code	CSB-PA973086
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NI60
Immunogen	Synthesized peptide derived from internal of Human ADCK3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ADCK3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Chaperone activity of bc1 complex-like; mitochondrial [Precursor]; Chaperone-ABC1-like; aarF domain-containing protein kinase 3; CABC1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COQ8A

Image



Western blot analysis of extracts from HepG2 cells and mouse muscle cells, using ADCK3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.