




LMO3 Antibody

Product Code	CSB-PA277951
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAP4
Immunogen	Synthesized peptide derived from internal of Human LMO3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total LMO3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	LIM domain only protein 3; Neuronal-specific transcription factor DAT1; Rhombotin-3; LMO3; RBTN3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LMO3

Image



Western blot analysis of extracts from COS-7 cells, using LMO3 antibody.