




ORAOV1 Antibody

Product Code	CSB-PA446256
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV07
Immunogen	Synthesized peptide derived from internal of Human ORAV1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ORAV1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Oral cancer overexpressed protein 1; Tumor amplified and overexpressed sequence 1; ORAOV1; TAOS1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ORAOV1

Image



Western blot analysis of extracts from 293 cells and HeLa cells, using ORAV1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.